

Product Datasheet

Human SPTAN1 ELISA Kit (orb1201330)

Description	Human Alpha-Fodrin (SPTAN1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q13813
Sensitivity	0.064ng/mL
Entrez	6709
Expiration Date	Please enquire.