

Product Datasheet

Human Smad7 ELISA Kit (orb1201337)

Description	Human Mothers Against Decapentaplegic Homolog 7 (Smad7) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O15105
Sensitivity	0.056ng/mL
Entrez	4092
Expiration Date	Please enquire.