

Product Datasheet

Human SIRPg ELISA Kit (orb1201342)

Description	Human Signal Regulatory Protein Gamma (SIRPg) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9P1W8
Sensitivity	0.057ng/mL
Entrez	55423
Expiration Date	Please enquire.