

Product Datasheet

Human SCCA2 ELISA Kit (orb1201348)

Description	Human Squamous Cell Carcinoma Antigen 2 (SCCA2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P48594
Sensitivity	0.064ng/mL
Entrez	6318
Expiration Date	Please enquire.