

Product Datasheet

Human PTHrP ELISA Kit (orb1201361)

Description	Human Parathyroid Hormone Related Protein (PTHrP) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P12272
Sensitivity	5.9pg/mL
Entrez	5744
Expiration Date	Please enquire.