

Product Datasheet

Human PKBg ELISA Kit (orb1201389)

Description	Human Protein Kinase B Gamma (PKBg) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9Y243
Sensitivity	0.56ng/mL
Entrez	10000
Expiration Date	Please enquire.