

## Product Datasheet

### Human PECAM1 ELISA Kit (orb1201404)

|                            |                                                                      |
|----------------------------|----------------------------------------------------------------------|
| <b>Description</b>         | Human Platelet/Endothelial Cell Adhesion Molecule (PECAM1) ELISA Kit |
| <b>Reactivity</b>          | Human                                                                |
| <b>Range</b>               | 31.2-2000pg/mL                                                       |
| <b>Tested Applications</b> | ELISA                                                                |
| <b>Note</b>                | For research use only                                                |
| <b>Application notes</b>   | sandwich                                                             |
| <b>Sample Types</b>        | serum, plasma, tissue homogenates or other biological fluids.        |
| <b>Uniprot ID</b>          | <b>P16284</b>                                                        |
| <b>Sensitivity</b>         | 13.6pg/mL                                                            |
| <b>Entrez</b>              | <b>5175</b>                                                          |
| <b>Expiration Date</b>     | Please enquire.                                                      |