

Product Datasheet

Human NOS1 ELISA Kit (orb1201441)

Description	Human Nitric Oxide Synthase 1, Neuronal (NOS1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P29475
Sensitivity	0.053ng/mL
Entrez	4842
Expiration Date	Please enquire.