

Product Datasheet

Human MYH1 ELISA Kit (orb1201450)

Description	Human Myosin Heavy Chain 1, Skeletal Muscle, Adult (MYH1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P12882
Sensitivity	0.059ng/mL
Entrez	4619
Expiration Date	Please enquire.