

Product Datasheet

Human MPP2 ELISA Kit (orb1201462)

Description	Human Membrane Protein, Palmitoylated 2 (MPP2) ELISA Kit
Reactivity	Human
Range	0.625-40ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q14168
Sensitivity	0.16ng/mL
Entrez	4355
Expiration Date	Please enquire.