

Product Datasheet

Human LAMA1 ELISA Kit (orb1201489)

Description	Human Laminin Subunit Alpha-1 (LAMA1) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P25391
Sensitivity	31pg/mL
Entrez	284217
Expiration Date	Please enquire.