

Product Datasheet

Human ICAM3 ELISA Kit (orb1201528)

Description	Human Intercellular Adhesion Molecule 3 (ICAM3) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P32942
Sensitivity	0.4ng/mL
Entrez	3385
Expiration Date	Please enquire.