

Product Datasheet

Human HSP90b1 ELISA Kit (orb1201535)

Description	Human Heat Shock Protein 90kDa Beta 1 (HSP90b1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P14625
Sensitivity	0.126ng/mL
Entrez	7184
Expiration Date	Please enquire.