

Product Datasheet

Human GnRH ELISA Kit (orb1201568)

Description	Human Gonadotropin Releasing Hormone (GnRH) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P01148
Sensitivity	0.064ng/mL
Entrez	2796
Expiration Date	Please enquire.