

Product Datasheet

Human GAL9 ELISA Kit (orb1201580)

Description	Human Galectin 9 (GAL9) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	000182
Sensitivity	34pg/mL
Entrez	3965
Expiration Date	Please enquire.