

Product Datasheet

Human GAL8 ELISA Kit (orb1201581)

Description	Human Galectin 8 (GAL8) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	000214
Sensitivity	0.116ng/mL
Entrez	3964
Expiration Date	Please enquire.