

## Product Datasheet

### Human G6PD ELISA Kit (orb1201590)

|                                |                                                                  |
|--------------------------------|------------------------------------------------------------------|
| <b>Description</b>             | Human Glucose 6 Phosphate<br>Dehydrogenase (G6PD) ELISA Kit      |
| <b>Reactivity</b>              | Human                                                            |
| <b>Range</b>                   | 6.25-400ng/mL                                                    |
| <b>Tested<br/>Applications</b> | ELISA                                                            |
| <b>Note</b>                    | For research use only                                            |
| <b>Application notes</b>       | sandwich                                                         |
| <b>Sample Types</b>            | serum, plasma, tissue homogenates or<br>other biological fluids. |
| <b>Uniprot ID</b>              | <b>P11413</b>                                                    |
| <b>Sensitivity</b>             | 2.35ng/mL                                                        |
| <b>Entrez</b>                  | <b>2539</b>                                                      |
| <b>Expiration Date</b>         | Please enquire.                                                  |