

Product Datasheet

Human CXCR4 ELISA Kit (orb1201632)

Description	Human Chemokine C-X-C-Motif Receptor 4 (CXCR4) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P61073
Sensitivity	0.060ng/mL
Entrez	7852
Expiration Date	Please enquire.