

Product Datasheet

Human CX3CL1 ELISA Kit (orb1201635)

Description	Human Chemokine C-X3-C-Motif Ligand 1 (CX3CL1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P78423
Sensitivity	0.053ng/mL
Entrez	6376
Expiration Date	Please enquire.