

Product Datasheet

Human CHEM ELISA Kit (orb1201657)

Description	Human Chemerin (CHEM) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q99969
Sensitivity	0.058ng/mL
Entrez	5919
Expiration Date	Please enquire.