

## Product Datasheet

### Human CENPE ELISA Kit (orb1201662)

|                            |                                                                 |
|----------------------------|-----------------------------------------------------------------|
| <b>Description</b>         | Human Centromere Protein E (CENPE)<br>ELISA Kit                 |
| <b>Reactivity</b>          | Human                                                           |
| <b>Range</b>               | 0.312-20ng/mL                                                   |
| <b>Tested Applications</b> | ELISA                                                           |
| <b>Note</b>                | For research use only                                           |
| <b>Application notes</b>   | sandwich                                                        |
| <b>Sample Types</b>        | tissue homogenates, cell lysates or<br>other biological fluids. |
| <b>Uniprot ID</b>          | <a href="#">Q02224</a>                                          |
| <b>Sensitivity</b>         | 0.124ng/mL                                                      |
| <b>Entrez</b>              | <a href="#">1062</a>                                            |
| <b>Expiration Date</b>     | Please enquire.                                                 |