

Product Datasheet

Human CCR5 ELISA Kit (orb1201673)

Description	Human Chemokine C-C-Motif Receptor 5 (CCR5) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P51681
Sensitivity	6.5pg/mL
Entrez	1234
Expiration Date	Please enquire.