

Product Datasheet

Human CCR4 ELISA Kit (orb1201674)

Description	Human Chemokine C-C-Motif Receptor 4 (CCR4) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P51679
Sensitivity	0.063ng/mL
Entrez	1233
Expiration Date	Please enquire.