

Product Datasheet

Human Bid ELISA Kit (orb1201704)

Description	Human BH3 Interacting Domain Death Agonist (Bid) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P55957
Sensitivity	0.115ng/mL
Entrez	637
Expiration Date	Please enquire.