

Product Datasheet

Human bACE1 ELISA Kit (orb1201707)

Description	Human Beta-Site APP Cleaving Enzyme 1 (bACE1) ELISA Kit
Reactivity	Human
Range	0.25-16ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, cerebrospinal fluid or other biological fluids.
Uniprot ID	P56817
Sensitivity	0.11ng/mL
Entrez	23621
Expiration Date	Please enquire.