

Product Datasheet

Human ANGPTL1 ELISA Kit (orb1201721)

Description	Human Angiopoietin Like Protein 1 (ANGPTL1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O95841
Sensitivity	0.055ng/mL
Entrez	9068
Expiration Date	Please enquire.