

Product Datasheet

Human ALT ELISA Kit (orb1201727)

Description	Human Alanine Aminotransferase (ALT) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P24298
Sensitivity	0.132ng/mL
Entrez	2875
Expiration Date	Please enquire.