

Product Datasheet

Human aHSP ELISA Kit (orb1201732)

Description	Human Alpha-Hemoglobin Stabilizing Protein (aHSP) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9NZD4
Sensitivity	0.055ng/mL
Entrez	51327
Expiration Date	Please enquire.