

Product Datasheet

Human AGER ELISA Kit (orb1201734)

Description	Human Advanced Glycosylation End Product Specific Receptor (AGER) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q15109
Sensitivity	0.058ng/mL
Entrez	177
Expiration Date	Please enquire.