

Product Datasheet

Human a1AGP ELISA Kit (orb1201749)

Description	Human Alpha-1-Acid Glycoprotein (a1AGP) ELISA Kit
Reactivity	Human
Range	3.12-200ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P02763
Sensitivity	1.27ng/mL
Entrez	5004
Expiration Date	Please enquire.