

Product Datasheet

Rat VEGFR1 ELISA Kit (orb1201753)

Description	Rat Vascular Endothelial Growth Factor Receptor 1 (VEGFR1) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P53767
Sensitivity	0.056ng/mL
Entrez	54251
Expiration Date	Please enquire.