

Product Datasheet

Rat Mlg ELISA Kit (orb1201760)

Description	Rat Monokine Induced By Interferon Gamma (Mlg) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q8K4B1
Sensitivity	0.057ng/mL
Entrez	246759
Expiration Date	Please enquire.