

Product Datasheet

Rat LFA1a ELISA Kit (orb1201762)

Description	Rat Lymphocyte Function Associated Antigen 1 Alpha (LFA1a) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q3T1L6
Sensitivity	0.063ng/mL
Entrez	308995
Expiration Date	Please enquire.