

Product Datasheet

Rat IRAK4 ELISA Kit (orb1201764)

Description	Rat Interleukin 1 Receptor Associated Kinase 4 (IRAK4) ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	D4A7K4
Sensitivity	5.9pg/mL
Entrez	300177
Expiration Date	Please enquire.