

Product Datasheet

Mouse TNFSF12 ELISA Kit (orb1201784)

Description	Mouse Tumor Necrosis Factor Ligand Superfamily, Member 12 (TNFSF12) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	O54907
Sensitivity	31.6pg/mL
Entrez	21944
Expiration Date	Please enquire.