

Product Datasheet

Guinea TGFb1 ELISA Kit (orb1201788)

Description	Guinea pig Transforming Growth Factor Beta 1 (TGFb1) ELISA Kit
Reactivity	Guinea pig
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, platelet-poor plasma, tissue homogenates, cell culture supernates or other biological fluids
Uniprot ID	Q9Z1Y6
Sensitivity	7.6pg/mL
Entrez	100135628
Expiration Date	Please enquire.