

Product Datasheet

Mouse PDGFB ELISA Kit (orb1201792)

Description	Mouse Platelet Derived Growth Factor Subunit B (PDGFB) ELISA Kit
Reactivity	Mouse
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P31240
Sensitivity	11.2pg/mL
Entrez	18591
Expiration Date	Please enquire.