

Product Datasheet

Canine IP10 ELISA Kit (orb1201795)

Description	Canine Interferon Gamma Induced Protein 10kDa (IP10) ELISA Kit
Reactivity	Canine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q5KSV9
Sensitivity	6.8pg/mL
Entrez	478432
Expiration Date	Please enquire.