

Product Datasheet

Monkey IL8 ELISA Kit (orb1201796)

Description	Monkey Interleukin 8 (IL8) ELISA Kit
Reactivity	Monkey
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P67813
Sensitivity	5.3pg/mL
Entrez	613028
Expiration Date	Please enquire.