

Product Datasheet

Guinea IL5 ELISA Kit (orb1201799)

Description	Guinea pig Interleukin 5 (IL5) ELISA Kit
Reactivity	Guinea pig
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	008987
Sensitivity	6.4pg/mL
Entrez	100135558
Expiration Date	Please enquire.