

Product Datasheet

Mouse IL34 ELISA Kit (orb1201804)

Description	Mouse Interleukin 34 (IL34) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q8R1R4
Sensitivity	5.3pg/mL
Entrez	76527
Expiration Date	Please enquire.