

Product Datasheet

Mouse IL24 ELISA Kit (orb1201806)

Description	Mouse Interleukin 24 (IL24) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q925S4
Sensitivity	6.2pg/mL
Entrez	93672
Expiration Date	Please enquire.