

Product Datasheet

Mouse IL22Ra2 ELISA Kit (orb1201807)

Description	Mouse Interleukin 22 Receptor Alpha 2 (IL22Ra2) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q80XF5
Sensitivity	0.119ng/mL
Entrez	237310
Expiration Date	Please enquire.