

Product Datasheet

Canine IL1a ELISA Kit (orb1201813)

Description	Canine Interleukin 1 Alpha (IL1a) ELISA Kit
Reactivity	Canine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O46612
Sensitivity	6.6pg/mL
Entrez	403782
Expiration Date	Please enquire.