

Product Datasheet

Canine IL17 ELISA Kit (orb1201817)

Description	Canine Interleukin 17 (IL17) ELISA Kit
Reactivity	Canine
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	C6L8D7
Sensitivity	31pg/mL
Entrez	481837
Expiration Date	Please enquire.