

Product Datasheet

Mouse FGF8 ELISA Kit (orb1201831)

Description	Mouse Fibroblast Growth Factor 8, Androgen Induced (FGF8) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P37237
Sensitivity	6.6pg/mL
Entrez	14179
Expiration Date	Please enquire.