

## Product Datasheet

### Mouse Bcl2L ELISA Kit (orb1201844)

|                            |                                                              |
|----------------------------|--------------------------------------------------------------|
| <b>Description</b>         | Mouse B-Cell CLL/Lymphoma 2 Like Protein (Bcl2L) ELISA Kit   |
| <b>Reactivity</b>          | Mouse                                                        |
| <b>Range</b>               | 78.1-5000pg/mL                                               |
| <b>Tested Applications</b> | ELISA                                                        |
| <b>Note</b>                | For research use only                                        |
| <b>Application notes</b>   | sandwich                                                     |
| <b>Sample Types</b>        | tissue homogenates, cell lysates or other biological fluids. |
| <b>Uniprot ID</b>          | <b>Q64373</b>                                                |
| <b>Sensitivity</b>         | 36.37pg/mL                                                   |
| <b>Entrez</b>              | <b>12048</b>                                                 |
| <b>Expiration Date</b>     | Please enquire.                                              |