

Product Datasheet

Rat SELP ELISA Kit (orb1201886)

Description	Rat Selectin, Platelet (SELP) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P98106
Sensitivity	0.138ng/mL
Entrez	25651
Expiration Date	Please enquire.