

Product Datasheet

Human SLAMF7 ELISA Kit (orb1201912)

Description	Human Signaling Lymphocytic Activation Molecule Family, Member 7 (SLAMF7) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9NQ25
Sensitivity	0.053ng/mL
Entrez	57823
Expiration Date	Please enquire.