

Product Datasheet

Human IL24 ELISA Kit (orb1201929)

Description	Human Interleukin 24 (IL24) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q13007
Sensitivity	27pg/mL
Entrez	11009
Expiration Date	Please enquire.